Van: [redacted], [redacted] <[redacted]@merck.com>
Verzonden: 11-02-2025 14:28
Aan: Mailbox ODBN <info@odbn.nl>
CC: [redacted] <[redacted]@odbn.nl>,
[redacted], [redacted] <[redacted]@aaenmaas.nl>
Onderwerp: Neomycine MSD Boxmeer

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Beste [redacted]

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[redacted]
[redacted]
MSD Animal Health

M. +316[redacted]

MSD Animal Health
The Science of Healthier Animals®

This e-mail message, together with any attachments, contains information of Merck & Co., Inc. (126 East Lincoln Ave., P.O. Box 2000, Rahway, NJ USA 07065) and/or its affiliates, that may be confidential, proprietary copyrighted and/or legally privileged. (Direct contact information for affiliates is available at - Contact us - MSD.) It is intended solely for the use of the individual or entity named on this message. If you are not the intended recipient, and have received this message in error, please notify us immediately by reply e-mail and then delete it from your system.

Method

Method:

Instrument parameter	Setting
Column	Thermo Hypersil Gold C18, 100x2.1 mm, 1.9 µm
Eluent A	0.1% HFBA* in water
Eluent B	0.1% HFBA* in methanol
Flow rate	0.3 mL/min
Column temperature	25 °C
Sample plate temperature	10 °C
Injection volume	100 µL
Runtime	11 min

*HFBA = heptafluorobutyric acid

Gradient:

Time (min)	% A	% B
0	90	10
6	5	95
8	5	95
8.5	90	10
11	90	10

QDa settings: positive mode


Ion	m/z	CV
[M+2H]2+	308.22	5
[M+H]+	615.41	15

Calibration line: 0-500 µg/L
SST/Limit solution: 10 µg/L

Sample preparation: To 4.5 mL sample add 0.5 mL 1% HFBA* in methanol.
For the spike samples: spike the solution with 50 µL Neomycin 1000 µg/L. Final concentration 10 µg/L.

Chromatogram of limit solution


Concentration of 10 µg/L
Ch1: m/z 308

S/N >> 10.
Method is suitable to detect levels at 10 µg/L in a clear solution without disturbing factors.

Chromatogram samples A0130 before evaporation

15-Apr-24 A0130
Ch1: m/z 308


Samples are measured and calculated according to the calibration line.

All samples show a peak for neomycin above >10 µg/L. Some samples show a peak above the calibration line (>1500 µg/L).

Chromatogram samples A0191 after evaporation

15-Apr-24 A0191
Ch1: m/z 308


Samples are not reproducible.

Chromatogram shows many disturbing peaks. This interferes the analysis.

Comparison with previous analysis: A0191 after evaporation

Ch1: m/z 308


Current chromatograms not comparable with previous chromatograms of samples.

There are a lot of interfering peaks present in current samples.

Chromatogram samples A0191 different channel m/z 615


Ch2: m/z 615

Chromatogram has less interfering peaks at the Rt of Neomycin than the channel m/z 308.

Channel m/z 615 does not match the results of channel m/z 308.


Limit determination based on m/z 615

- Less disturbing factors around Neomycine at m/z 615
- The peak of Neomyice in all samples is much lower than the 10 ug/L peak in the standard (black line) as can be seen in the chromatogram.
- From these chromatograms it can be concluded that the amount of neomycine in the samples is below 10 ug/L based on the results obtained from m/z 615

Limit determination based on m/z 615

Samples A0191	[M+H]+ (m/z 615)	
	Spike Rec. (%)	Result (µg/L)
Blank	102	<10
Apr	82	<10
May	89	<10
Jun	69	<10
Jul	63	<10
Aug	71	<10
Sep	88	<10

- In some of the samples a low recovery is found
- Nevertheless, the neomycine peak in the samples is much lower than in spiked samples (10 ug/L)
- This indicates that the amount neomycine in the samples is much lower than 10 ug/L

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